




Create a Service Code to unlock your Valve Tester V2 device.

1. Download and install firmware.

https://opensys.com.ua/system/storage/download/kYXm-ValveTester_v1_97_20250310_2.0_11.3_NL0_DEF_LIM_ENABLE_0.exe


2. Connect the device to the AC (110/260 V 50/60 Hz) network and to the computer via the USB cable. Turn on the device.
3. Update firmware Valve Tester



4. Download and install software for Valve Tester.


<https://opensys.com.ua/system/storage/download/46Eu-ValveTester-1.3.0.msi>

Start software Valve Tester






6. Connect the device using the Valvetester program



7. Create Service Request



8. Contact your seller.
9. After generating the service code, your seller will send you a file with the code. It will be necessary to press the "execute service code" button and select the file with the code you received

